

TABS® Basic

Model 25000



Features

- ▲ **Limited Lifetime Warranty***
- ▲ **Easy to use and install**
- ▲ **Durable, ergonomic design**
- ▲ **Easy transfer between bed, wheelchair and room**
- ▲ **Two wall mounting choices**
- ▲ **Audible low battery indication**
- ▲ **Locking garment clip**
- ▲ **No Velcro®**

Options

- ▲ **Bed bracket**
- ▲ **Over-the-door bracket**
- ▲ **Various pull cord lengths**

Specifications

- ▲ **Weight: Not more than 2 lbs.**
- ▲ **Sound pressure: >75dBA <90dBA**
- ▲ **Power: 9V alkaline battery**
- ▲ **Shunt pull force: Min. 10 oz.; max. 3 lbs.**

The Basic Tool for Fall Management

The TABS Basic mobility monitor is a restraint-free medical device that:

- **Alerts you when a patient is leaving a bed, chair, wheelchair or room**
- **May serve as a reminder to patients to call for assistance before rising**
- **Helps to reduce the risk of falls**
- **Helps improve staff response time**

With a simple disk, clip, and cord system, the TABS Basic monitor allows you to determine a personal safety zone for each patient and then alerts you when patients leave their personal safety zones.

TABS pressure pads are not compatible with the TABS Basic model. Other TABS monitors are available that can use pressure pads. Call your TABS distributor to determine which monitor is best suited for your purpose.

From the rugged case and secure mounting bracket to the locking garment clip, TABS products are among the highest quality fall prevention monitors available. Their versatility allows you to easily transfer a single monitor between bed, chair, wheelchair, toilet, and door. The over-the-door bracket enables use as a room monitor.



TABS monitors comply with JCAHO Patient Safety Standards. There are no on/off, volume control or delay switches, which helps to minimize human error.

Distributed by: