

TABS® Select

Model 25011



Features

- ▲ **Limited Lifetime Warranty***
- ▲ **String and/or pressure pad use**
- ▲ **Four distinct alarm tones**
- ▲ **Resets during pressure pad use when patient weight returns to pad**
- ▲ **Audible low battery indication**
- ▲ **Two wall mounting choices**
- ▲ **Locking garment clip**
- ▲ **No Velcro®**

Options

- ▲ **Nurse call interface**
- ▲ **Bed bracket**
- ▲ **See Accessories page for more**

Specifications

- ▲ **Weight: Not more than 2 lbs.**
- ▲ **Sound pressure: >75dBA <90dBA**
- ▲ **Power: 9V alkaline battery**
- ▲ **Shunt pull force: Min. 10 oz.; max. 3 lbs.**

The Select Tool for Fall Management

The TABS Select monitor is a restraint-free medical device that allows you to:

- **Choose the appropriate monitoring methods for each individual. Use the pull-string assembly, pressure pad, or both simultaneously.**
- **Know when a patient is leaving a bed, chair, wheelchair, or toilet.**
- **Identify a single patient or an entire wing with separate tones. Four distinct alarm tones are available.**
- **Avoid disturbances. With the optional nurse call interface, you can set the alarm to sound only at the nurse station and remain silent in the room.**

From the rugged case and secure mounting bracket to the locking garment clip, TABS products are among the highest quality fall prevention monitors available. Their versatility allows you to easily transfer a single monitor between bed, chair, wheelchair, and toilet.



TABS monitors comply with JCAHO Patient Safety Standards. There are no on/off, volume control or delay switches, which helps to minimize human error.

Distributed by:

800-824-4490 | Fax 402-475-4023 | www.stanleyhealthcare.com

STANLEY
Healthcare Solutions